

Product Datasheet

Fluorescent Mid Range DNA Ladder (orb1734251)

Description	Fluorescent Mid Range DNA Ladder consists of 12 DNA fragments ranging from 100 bp to 3 kb. It is des...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Tartrazine and Xylene Cyanol FF; Color: green
Concentration	130 ng/μl
Storage	store at -20 °C. store dark. Short term storage (up to 1 month) at 4 °C possible
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.