

Product Datasheet

1 kb DNA Ladder (orb1734253)

Description	1 kb DNA Ladder consist of 10 DNA fragments ranging from 1 kb to 10 kb and is supplied in ready-to-u...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, Bromophenol Blue and Xylene Cyanol FF; Color: blue
Concentration	105 ng/μl
Storage	store at -20 °C. Short term storage (up to 1 month) at 4 °C possible
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.