

Product Datasheet

200 bp DNA Ladder (orb1734255)

Description	200 bp DNA Ladder consist of 10 DNA fragments ranging from 200 bp to 2 kb and is supplied in ready-t...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, Bromophenol Blue and Xylene Cyanol FF; Color: blue
Concentration	105 ng/μl
Storage	store at -20 °C. Short term storage (up to 1 month) at 4 °C possible
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.