

## Product Datasheet

### Low Range DNA Ladder (orb1734262)

<b>Description</b>	Low Range DNA Ladder consists of 11 DNA fragments ranging from 50 bp to 1 kb. It is designed to show...
<b>Form/Appearance</b>	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 8% (w/v) Glycerol, Tartrazine and Xylene Cyanol FF; Color: green
<b>Concentration</b>	120 ng/μl
<b>Storage</b>	store at -20 °C. Short term storage (up to 1 month) at 4 °C possible
<b>Note</b>	For research use only
<b>Application notes</b>	<b>Recommended Load:</b> 5 μl per lane.
<b>Expiration Date</b>	12 months from date of receipt.