

## Product Datasheet

### Azide-PEG3-FLAG (orb1734977)

<b>Description</b>	Azide-PEG 3-FLAG enables the FLAG-tag (DYKDDDDK) attachment to any terminal Alkyne-or Cyclooctyne (e...
<b>Form/Appearance</b>	solid; Color: white to off-white
<b>Storage</b>	store at -20 °C. Short term exposure (up to 1 week cumulative) to ambient temperature possible
<b>Note</b>	For research use only
<b>Formula</b>	$C_{59}H_{88}N_{16}O_{27}S$
<b>Purity</b>	≥ 95% (HPLC)
<b>MW</b>	Theoretical MW: 1485.49 g/mol; Detected MW: 1484.57 g/mol
<b>Solubility (25°C)</b>	water, PBS (up to 30 mM tested)
<b>Expiration Date</b>	12 months from date of receipt.