

Product Datasheet

Tris-HCl Buffer 1 M pH 8.5 (orb1734979)

Description	Tris-HCl Buffer is included in PCR reaction buffer, used for in vitro cell culture, enzyme assays an...
Form/Appearance	liquid; pH: 8.5 ± 0.1 (22 °C)
Concentration	1 M
Storage	store at 4 °C
Note	For research use only
Formula	$C_4H_{11}NO_3 \cdot HCl$
MW	Theoretical MW: 157.60 g/mol
CAS Number	1185-53-1
Expiration Date	12 months from date of receipt.