



Product Datasheet

STAT3 Phospho (Tyr705) Antibody (orb1735066)

Description	STAT3 Phospho (Tyr705) Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to residues surrounding Tyr705 of human STAT3
Target	STAT3 Phospho (Tyr705), monoclonal
Storage	storage is recommended at 4°C. Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200-1:400
Expiration Date	12 months from date of receipt.



orb1735066
Anti-
Phospho
STAT3
(Tyr705)
m...



Serum
starved
A431 cells
were
untreated ...