

## Product Datasheet

### Norovirus NS6 Antibody (orb1735086)

<b>Description</b>	Norovirus NS6 Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, WB
<b>Target</b>	Norovirus NS4 antibody
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	Western Blot 1:3000-1:15000 Immunofluorescence 1:500 - 1:1000 ELISA (1:5000-10000)
<b>Expiration Date</b>	12 months from date of receipt.

