

Product Datasheet

Norovirus VP1 Antibody (orb1735088)

Description	Norovirus VP1 Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	Norovirus VP1 Monoclonal antibody
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	Western Blot 1:3000-1:15000 Immunofluorescence 1:500 - 1:1000 ELISA (1:5000-10000)
Expiration Date	12 months from date of receipt.

