

## Product Datasheet

# Anti-HRSV-A/B F Protein Antibody (ADI-15600) (orb1735115)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Anti-HRSV-A/B F Protein Antibody<br>(ADI-15600)  |
| <b>Species/Host</b>    | Mouse  |
| <b>Reactivity</b>      | Virus  |
| <b>Conjugation</b>     | Unconjugated   |
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up<br>to 2 weeks. For long term storage<br>store at -20°C in small aliquots to<br>prevent freeze-thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Isotype</b>         | IgG2a  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |