

Product Datasheet

Anti-HPV16 L2/Minor capsid protein L2 Antibody (MAb24B) (orb1735117)

Description	Anti-HPV16 L2/Minor capsid protein L2 Antibody (MAb24B)
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1, kappa
Expiration Date	12 months from date of receipt.