

## Product Datasheet

# Anti-Human parvovirus B19/HPV-B19 VP2/Viral protein 2 Antibody (SAA0341) (orb1735118)

<b>Description</b>	Anti-Human parvovirus B19/HPV-B19 VP2/Viral protein 2 Antibody (SAA0341)
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1, kappa
<b>Expiration Date</b>	12 months from date of receipt.