

Product Datasheet

Anti-HPV16 E6/Protein E6 Antibody (Nb9) (orb1735119)

Description	Anti-HPV16 E6/Protein E6 Antibody (Nb9)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	VHH-mFC
Expiration Date	12 months from date of receipt.