

Product Datasheet

Anti-HPV16 E7/Protein E7 Antibody (Nb4) (orb1735121)

Description	Anti-HPV16 E7/Protein E7 Antibody (Nb4)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	VHH-mFC
Expiration Date	12 months from date of receipt.