

Product Datasheet

Anti-HPV16 L1/Major capsid protein L1 Antibody (H16.V5) (orb1735127)

Description	Anti-HPV16 L1/Major capsid protein L1 Antibody (H16.V5)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Preservatives	PBS buffer, pH7.4
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C for 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1, kappa
Uniprot ID	P03101
Expiration Date	12 months from date of receipt.