

Product Datasheet

Anti-HRSV-A/B G/Major surface glycoprotein G Antibody (131-2G) (orb1735132)

Description	Anti-HRSV-A/B G/Major surface glycoprotein G Antibody (131-2G)
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1, kappa
Expiration Date	12 months from date of receipt.