

Product Datasheet

Rat IgG2a kappa Anti-Mouse IgM heavy chain (HRP) 1 mg (orb1735446)

Description	Rat IgG2a kappa Anti-Mouse IgM heavy chain (HRP) 1 mg
Species/Host	Rat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	CS, ELISA, ELISpot, IA, IAC
Target	Mouse IgM heavy chain
Preservatives	PBS/50%glycerol. PBS/50%glycerol
Form/Appearance	HRP
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2a kappa
Clone Number	LO-MM-9
Antibody Type	Secondary Antibody
Uniprot ID	P01872
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)