

Product Datasheet

Rat IgM kappa Anti-Mouse IgG3 heavy chain (FITC) 0,5 mg (orb1735509)

Description	Rat IgM kappa Anti-Mouse IgG3 heavy chain (FITC) 0,5 mg
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	ELISA, IA
Target	Mouse IgG3 heavy chain
Preservatives	PBS/0,1%NaN ₃ /50%glycérol. PBS/0,1%NaN ₃ /50%glycérol
Form/Appearance	FITC
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgM kappa
Clone Number	LO-MG3-7
Antibody Type	Secondary Antibody
Uniprot ID	P03987
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)