

Product Datasheet

His Tag (N-Ter & C-Ter) Antibody (HRP) (orb1735759)

Description	His Tag (N-Ter & C-Ter) Antibody (HRP)
Species/Host	Rat
Conjugation	HRP
Tested Applications	ELISA, WB
Target	His Tag (N-Ter & C-Ter)
Preservatives	PBS/50%glycerol. PBS/50%glycerol
Form/Appearance	HRP
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2b kappa
Expiration Date	12 months from date of receipt.