

Product Datasheet

Rat IgG2a kappa (Isotype control) Antibody (orb1736020)

Description	Rat IgG2a kappa (Isotype control) Antibody
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IA
Preservatives	PBS/0,1%NaN3. PBS/0,1%NaN3
Form/Appearance	Purified
Concentration	1mg/mL
Storage	-20°C
Note	For research use only
Isotype	IgG2a kappa
Expiration Date	12 months from date of receipt.