

## Product Datasheet

# Rat IgG2a kappa (Isotype control) Antibody (orb1736023)

<b>Description</b>	Rat IgG2a kappa (Isotype control) Antibody
<b>Species/Host</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IA
<b>Preservatives</b>	PBS/0,1%NaN3. PBS/0,1%NaN3
<b>Form/Appearance</b>	Purified
<b>Concentration</b>	1mg/mL
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a kappa
<b>Expiration Date</b>	12 months from date of receipt.