

Product Datasheet

Rat IgA kappa (Isotype control) Antibody (orb1736039)

| | |
|----------------------------|--|
| Description | Rat IgA kappa (Isotype control) Antibody |
| Species/Host | Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IA |
| Preservatives | PBS/0,1%NaN3. PBS/0,1%NaN3 |
| Form/Appearance | Purified |
| Concentration | 1mg/mL |
| Storage | -20°C |
| Note | For research use only |
| Isotype | IgA kappa |
| Expiration Date | 12 months from date of receipt. |