

Product Datasheet

anti-VEGFR2 VHH antibody (orb1736364)

| | |
|----------------------------|--|
| Description | Primary VHH single domain antibody |
| Species/Host | Camelus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Target | Human Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) |
| Preservatives | 25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin |
| Form/Appearance | Liquid |
| Concentration | 1mg/mL |
| Storage | Storage Condition: 4°C; Do not freeze. Shelf Life: 1 year |
| Note | For research use only |
| Isotype | Camelid VHH |
| Clonality | Monoclonal |
| Purity | > 95% as determined by SDS-PAGE |
| MW | 15 kDa |
| Expiration Date | 12 months from date of receipt. |