

## Product Datasheet

# Human Glial Fibrillary Acidic Protein (GFAP) Antibody (orb1736407)

<b>Description</b>	Primary VHH single domain antibody
<b>Species/Host</b>	Camelus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Target</b>	Human Glial Fibrillary Acidic Protein (GFAP)
<b>Preservatives</b>	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Storage Condition: 4°C; Do not freeze. Shelf Life: 1 year
<b>Note</b>	For research use only
<b>Isotype</b>	Camelid VHH
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95% as determined by SDS-PAGE
<b>MW</b>	15 kDa
<b>Expiration Date</b>	12 months from date of receipt.