

Product Datasheet

Human Tumor Necrosis Factor Receptor 1 (hTNFR1) Antibody (orb1736426)

Description	Primary VHH single domain antibody
Species/Host	Camelus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Target	Human Tumor Necrosis Factor Receptor 1 (hTNFR1)
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin
Form/Appearance	Liquid
Concentration	1mg/mL
Storage	Storage Condition: 4°C; Do not freeze. Shelf Life: 1 year
Note	For research use only
Isotype	Camelid VHH
Clonality	Monoclonal
Purity	> 95% as determined by SDS-PAGE
MW	15 kDa
Expiration Date	12 months from date of receipt.