

## Product Datasheet

# Pig VEGFD(Vascular Endothelial Growth Factor D) ELISA Kit (orb1736459)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Porcine
<b>Range</b>	31.25-2000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	11.8 pg/mL
<b>Expiration Date</b>	Please enquire.

