

## Product Datasheet

### Mouse TRKB (Tyrosine Kinase B) ELISA Kit (orb1736470)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Mouse
<b>Range</b>	93.75-6000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Tissue homogenates, cell lysates and other biological fluids.
<b>Sensitivity</b>	37.24 pg/mL
<b>Expiration Date</b>	Please enquire.

