

## Product Datasheet

# Chicken VDR(Vitamin D Receptor) ELISA Kit (orb1736480)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Gallus
<b>Range</b>	0.63-40 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids
<b>Uniprot ID</b>	<a href="#">O42392</a>
<b>Sensitivity</b>	0.25 ng/mL
<b>Expiration Date</b>	Please enquire.

