

## Product Datasheet

# Mouse CD200R1(CD200 Receptor 1) ELISA Kit (orb1736515)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Mouse
<b>Range</b>	78.13-5000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates and other biological fluids
<b>Uniprot ID</b>	<a href="#">Q9ES57</a>
<b>Sensitivity</b>	27.2 pg/mL
<b>Expiration Date</b>	Please enquire.

