

Product Datasheet

Mouse TYROBP(TYRO Protein Tyrosine Kinase Binding Protein) ELISA Kit (orb1736517)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.32-20 ng/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	054885
Sensitivity	0.16 ng/mL
Expiration Date	Please enquire.

