

## Product Datasheet

### Goat T3(Triiodothyronine) ELISA Kit (orb1736585)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Goat
<b>Range</b>	0.16-10 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma and other biological fluids
<b>Sensitivity</b>	0.047 ng/mL
<b>Expiration Date</b>	Please enquire.

