

Product Datasheet

Cattle CXCL6(Chemokine C-X-C-Motif Ligand 6) ELISA Kit (orb1736590)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Bovine
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P80221
Sensitivity	6.47 pg/mL
Expiration Date	Please enquire.

