

## Product Datasheet

# Human BUB1B(Mitotic Checkpoint Serine/threonine-Protein Kinase BUB1 Beta) ELISA Kit (orb1736598)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	0.32-20 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma and other biological fluids
<b>Uniprot ID</b>	<b>O60566</b>
<b>Sensitivity</b>	0.12 ng/mL
<b>Expiration Date</b>	Please enquire.

