

## Product Datasheet

### Rat MT(Melatonin) ELISA Kit (orb1736603)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Rat
<b>Range</b>	7.82-500 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids
<b>Sensitivity</b>	2.81 pg/mL
<b>Expiration Date</b>	Please enquire.

