

## Product Datasheet

# Rat PDL1(Programmed Cell Death Protein 1 Ligand 1) ELISA Kit (orb1736657)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Rat
<b>Range</b>	78.13-5000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids
<b>Sensitivity</b>	31.2 pg/mL
<b>Expiration Date</b>	Please enquire.

