

## Product Datasheet

# Human NEFH(Neurofilament Heavy Polypeptide) ELISA Kit (orb1736661)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	15.63-1000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids
<b>Uniprot ID</b>	<b>P12036</b>
<b>Sensitivity</b>	6.81 pg/mL
<b>Expiration Date</b>	Please enquire.

