

Product Datasheet

Chicken CCK(Cholecystikin) ELISA Kit (orb1736686)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Gallus
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9PU41
Sensitivity	4.82 pg/mL
Expiration Date	Please enquire.

