

Product Datasheet

Human Siglec5(Sialic acid-binding Ig-like lectin 5) ELISA Kit (orb1736691)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	78.13-5000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O15389
Sensitivity	27.8 pg/mL
Expiration Date	Please enquire.

