

Product Datasheet

Human Siglec9(Sialic acid-binding Ig-like lectin 9) ELISA Kit (orb1736692)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.47-30 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9Y336
Sensitivity	0.18 ng/mL
Expiration Date	Please enquire.

