

## Product Datasheet

# Mouse SLC6A1(Sodium- and chloride-dependent GABA transporter 1) ELISA Kit (orb1736722)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 0.16-10 ng/mL

**Note** For research use only

**Sample Types** serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

**Uniprot ID** [P31648](#)

**Sensitivity** 0.074 ng/mL

**Expiration Date** Please enquire.

