

## Product Datasheet

# Human C1QTNF7(C1q And Tumor Necrosis Factor Related Protein 7) ELISA Kit (orb1736739)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	2.5-160 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9BXJ2</a>
<b>Sensitivity</b>	1.36 ng/mL
<b>Expiration Date</b>	Please enquire.

