

## Product Datasheet

# Guinea pig C3(Complement Component 3) ELISA Kit (orb1736740)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Guinea pig
<b>Range</b>	0.32-20 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids.
<b>Sensitivity</b>	0.14 ng/mL
<b>Expiration Date</b>	Please enquire.

