

Product Datasheet

Cat CXCR4(Chemokine C-X-C-Motif Receptor 4) ELISA Kit (orb1736771)

| | |
|------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Feline |
| Range | 78.13-5000 pg/mL |
| Note | For research use only |
| Sample Types | tissue homogenates, cell lysates and other biological fluids |
| Uniprot ID | P56498 |
| Sensitivity | 36.8 pg/mL |
| Expiration Date | Please enquire. |

