

Product Datasheet

Human Recombinant ACE-2 (Fc tag) (orb1736820)

Description	Recombinant human ACE-2 tag fusion fused to IgG1 Fc. Produced in Nicotiana benthamiana plants and fo...
Species/Host	Plant
Tested Applications	ELISA, WB
Preservatives	0.1 M Phosphate Buffered Saline (PBS), pH = 7.4
Form/Appearance	Liquid
Concentration	1,0 mg/ml
Storage	-20°C Freezer
Note	For research use only
Formula	Monomer
Isotype	IgG1
Clonality	Monoclonal
Purity	92.03 % as determined by Mass Spectrometry.
Source	Human
MW	± 220 kDa
Uniprot ID	Q9BYF1
Hazard Information	Non-Hazardous
Dilution Range	Optimal dilutions/concentrations should be determined by the user.
Expiration Date	6 months from date of receipt.