

---

## Product Datasheet

### Native Human hCG protein Antibody (orb1736828)

**Description**

Recombinant rabbit monoclonal antibody

protein. This product i...

**Species/Host**

Plant

**Conjugation**

Unconjugated

**Tested**

ELISA

**Applications**
**Target**

Native Human hCG protein

**Preservatives**

 0.1 M Phosphate Buffered Saline (PBS),  
pH = 7.6

**Form/Appearance**

Liquid

**Concentration**

1,0 mg/ml

**Storage**

-20°C Freezer

**Note**

For research use only

**Formula**

Full-length IgG

**Isotype**

IgG1

**Clonality**

Monoclonal

**Purity**

≥ 90 % as determined by SDS-PAGE

**Source**

Rabbit

**MW**

± 150 kDa

**Uniprot ID**
**PODN86**
**Sensitivity**

 Refer to ELISA dose response graph in  
Datasheet for ELISA sensitivity.

**Hazard**

Non-Hazardous

**Information**
**Dilution Range**

 Suggested dilution is 1: 3 000-1: 256 000  
for ELISA. Optimal  
dilutions/concentrations should be  
determined by the user.

**Expiration Date**

12 months from date of receipt.