

Product Datasheet

HIV-1 p24 Antibody (orb1736840)

Description

Recombinant human monoclonal antibody

is produced in ...

Species/Host

Plant

Conjugation

Unconjugated

Tested Applications

ELISA

Target

HIV-1 p24

Preservatives

0.1 M Phosphate Buffered Saline (PBS), pH = 7.4

Form/Appearance

Liquid

Concentration

1,0 mg/ml

Storage

-20°C Freezer

Note

For research use only

Formula

Full-length IgG

Isotype

IgG1

Clonality

Monoclonal

Purity

≥ 90 % as determined by SDS-PAGE.

Source

Human

MW

± 150 kDa

Uniprot ID**Q9WMV8****Sensitivity**

Refer to ELISA dose response graph in Datasheet for ELISA sensitivity.

Hazard Information

Non-Hazardous

Dilution Range

Suggested dilution is 1: 1 000 - 1: 10 000 for ELISA. Optimal dilutions/concentrations should be determined by the user.

Expiration Date

12 months from date of receipt.