

---

## Product Datasheet

**SARS-CoV-2 Spike protein (RBD of S1 subunit)  
Antibody (orb1736854)**

**Description**

Recombinant human monoclonal

	Binding Domain of SARS-CoV-2 Spike protei...
<b>Species/Host</b>	Plant
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	SARS-CoV-2 Spike protein (RBD of S1 subunit)
<b>Preservatives</b>	0.1 M Phosphate Buffered Saline (PBS), pH = 7.4
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	-20°C Freezer
<b>Note</b>	For research use only
<b>Formula</b>	Full-length IgG
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	99.30 % as determined by Mass Spectrometry.
<b>Source</b>	Human
<b>MW</b>	± 150 kDa
<b>Uniprot ID</b>	<a href="#">P59594</a>
<b>Hazard Information</b>	Non-Hazardous
<b>Dilution Range</b>	Optimal dilutions/concentrations should be determined by the user.
<b>Expiration Date</b>	12 months from date of receipt.