

---

## Product Datasheet

### **SARS-CoV-2 Spike protein (S1 Subunit) Antibody (orb1736855)**

**Description**

Recombinant Alpaca monoclonal

CoV-2 Spike protein. This product ...

**Species/Host**

Plant

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, WB

**Target**

SARS-CoV-2 Spike protein (S1 Subunit)

**Preservatives**

0.1 M Phosphate Buffered Saline (PBS), pH = 7.4

**Form/Appearance**

Liquid

**Concentration**

1.0 mg/ml

**Storage**

-20°C Freezer

**Note**

For research use only

**Formula**

Full-length IgG

**Isotype**

IgG2b

**Clonality**

Monoclonal

**Purity**

&gt; 85 % as determined by SDS-PAGE.

**Source**

Alpaca

**MW**

± 150 kDa

**Uniprot ID**
**P59594**
**Hazard Information**

Non-Hazardous

**Dilution Range**

Optimal dilutions/concentrations should be determined by the user.

**Expiration Date**

12 months from date of receipt.