

---

## Product Datasheet

### Human Interferon Gamma (IFN- $\gamma$ ) Antibody (orb1736856)

**Description**

Recombinant rabbit monoclonal

This product is produced in Nicotia...

**Species/Host**

Plant

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

WB

**Target**

 Human Interferon Gamma (IFN- $\gamma$ )

**Preservatives**

 0.1 M Phosphate Buffered Saline (PBS),  
pH = 7.4

**Form/Appearance**

Liquid

**Concentration**

1.0 mg/ml

**Storage**

-20°C Freezer

**Note**

For research use only

**Formula**

Full-length IgG

**Isotype**

IgG1

**Clonality**

Monoclonal

**Purity**

&gt; 90 % as determined by SDS-PAGE.

**Source**

Rabbit

**MW**
 $\pm$  150 kDa

**Uniprot ID**
**P01579**
**Sensitivity**

 Detected from as low as 31.25 ng for  
WB.

**Hazard Information**

Non-Hazardous

**Dilution Range**

 Suggested dilution is 1: 1 000 - 1: 7 500  
for WB. Optimal  
dilutions/concentrations should be  
determined by the user.

**Expiration Date**

12 months from date of receipt.