

---

## Product Datasheet

### Human Tubulin beta chain Antibody (orb1736858)

**Description**

Recombinant rabbit monoclonal antibody

produced in Nicotiana be...

**Species/Host**

Plant

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested Applications**

WB

**Target**

Human Tubulin beta chain

**Preservatives**

 0.1 M Phosphate Buffered Saline (PBS),  
pH = 7.4

**Form/Appearance**

Liquid

**Concentration**

1.0 mg/ml

**Storage**

-20°C Freezer

**Note**

For research use only

**Formula**

Full-length IgG

**Isotype**

IgG1

**Clonality**

Monoclonal

**Purity**

 97.67 % as determined by Mass  
Spectrometry.

**Source**

Rabbit

**MW**

± 150 kDa

**Uniprot ID**
**P07437**
**Sensitivity**

Detected from as low as 62.5 ng for WB.

**Hazard Information**

Non-Hazardous

**Dilution Range**

 Suggested dilution is 1: 1 000-1: 5 000  
for WB. Optimal dilutions/concentrations  
should be determined by the user.

**Expiration Date**

12 months from date of receipt.